


An Overview of Capabilities at the NIST Center for Neutron Research

Juscelino B. Leão
7th International Workshop on
Sample Environment at Neutron Facilities
Sydney, Australia

- The Tool Belt
- The Right Tool for the Right Job
- News
- Planning
- Safety vs. Results
- R&D

The Tool Belt


What Range?



<http://www.ncnr.nist.gov/equipment/ancequip.html>

The Right Tool for the Right Job

What is important? At what instrument?




Data Acquisition Time

Accuracy and Precision

Simultaneous Multi-Techniques

<http://www.ncnr.nist.gov/equipment/ancequip.html>

At what cost?



News



wide

Fridge

pendak

11.5 T we

Newly commissioned SANS


nty of v

sys

ma

gn

nt



y

g w



Planning

P_{rior}

P_{reparation}

P_{revents}

P_{oor}

P_{erformance}

- Preparation (S.E. consultation)
- Submission
- Review
 - ✓ Scientific – national & international experts
 - ✓ Technical & Safety Review (S.E.)
 - ✓ Proposal Assessment Committee (PAC)
- Scheduling
 - ✓ User Office
 - ✓ Instrument Scientists (S.E. consultation)
(24 Hrs. prior)
- Arrival
- Completion
- Customer feedback



<http://www.ncnr.nist.gov/>

Reality

- Preparation ~~(S.E. consultation)~~
- Submission
- Review
 - ✓ Scientific – national & international experts
 - ✓ Technical & Safety Review (BYO S.E or n/a)
 - ✓ PAC (Proposal Calls for Ambient Conditions)
- Scheduling **SE GROUP IS RELIEVED**
 - ✓ User Office
 - ✓ Instrument Scientists
- Arrival
 - ~~(S.E. consultation)~~
 - (Day of experiment)
 - (Parameters changed)
 - SE scrambles to meet demand
- Completion
- Customer feedback **@!%!**#~@#!! SE Group**


<http://www.ncnr.nist.gov/>

Safety vs. Results

Stick to the plan



Slow down and be vigilant



In doubt ask




Know your sample

Yes! My sample
IS dry!



^3He Polarization


^3He program for scattering applications
SEOP lab



Currently available for 3-axis,
reflectometry, and SANS

≈ 20 experiments/year

Computer Controlled Gas Handling Manifold



Pressures up to 200 bar

Resolution: 0.01% of F.S.

Accuracy: $\pm 1\%$ Reading


Flow restrictors to prevent gauge saturation, vacuum pump

Intrinsically safe wiring

Explosion proof enclosure, gauges and valves

Sample volume protected via expansion volume

Easy pull-down menus
Supports scripting for remote beam line experiments
Works either as a stand-alone program or in the
background controlled by ICP
Expandable



Software Screen Shot




Multi-stage SAmple Changer


Prototype:

6 Samples

$253\text{ K} \leq T \leq 393\text{ K}$


$t_{\Delta T=100\text{K}} : 1\text{ min. (0.25 K acc.)}$
 $: 5\text{ min. (0.1 K acc.)}$






Receives pre-existing SANS
demountable cells
Ease of sample change

Spring-loaded
plunger






Torlon vrs PTFE:

Thermal Conductivity 0.250 W/m-K

Deflection Temperature at 1.8 Mpa

PTFE: 93.3 °C

Torlon: 278 °C



SANS Rheometers



MCR 501

Rheology

Couette Geometry

Cone-Plate/Plate-Plate

Torque Ranges

MCR501: 0.1 μNm – 230 mNm

MCR301: 0.1 μNm – 200 mNm

Shear Stress: 0.5 mPa – 5.5 MPa



RheoSANS

Couette Geometry

1,3 and 2,3 planes

Cups and Bobs from Titanium and
Quartz

Static SANS Cell

Time Resolved Measurements



Sample Cans



<http://www.ncnr.nist.gov/equipment/ancequip.html>



Air sensitive/gas loading

$P_{\max} = 5 \text{ bar (V)}$

$4 \text{ K} \leq T \leq 800 \text{ K}$

Heated gas line available for
methane and CO_2

SANS Gas Adsorption

$P_{\max} = 1,000 \text{ bar}$

$\text{LN}_2 < T < 350 \text{ K}$

Beam divergence angle $\theta \approx 20^\circ$

